

## PREVENTIVE MAINTENANCE PROGRAM



## SAFE CleanBox STK 113 – STK 103

## PREVENTIVE MAINTENANCE PROGRAM

	1st YEAR		
	1st quarter	2nd quarter	3rd quarter
N. technicians	1	1	1
N. h/ps	2	2	3
INTERVENTION	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE
EXPECTED REPLACEMENTS			WIPER BRUSH ASPIRATOR + FILTERS

	2nd YEAR		
	1st quarter	2nd quarter	3rd quarter
N. technicians	1	1	1
N. h/ps	2	2	4
INTERVENTION	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE
EXPECTED REPLACEMENTS			WIPER BRUSH ASPIRATOR + FILTERS FLANGE GASKETS DOOR GASKETS

	3rd YEAR		
	1st quarter	2nd quarter	3rd quarter
N. technicians	1	1	1
N. h/ps	2	2	4
INTERVENTION	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE
EXPECTED REPLACEMENTS			WIPER BRUSH ASPIRATOR + FILTERS AIR/WATER TUBES- HANDPIECES

	4th YEAR		
	1st quarter	2nd quarter	3rd quarter
N. technicians	1	1	1
N. h/ps	2	2	4
INTERVENTION	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE
EXPECTED REPLACEMENTS			WIPER BRUSH EXHAUSTER + FILTERS FLANGE GASKETS DOOR GASKETS

	5th YEAR		
	1st quarter	2nd quarter	3rd quarter
N. technicians	1	1	1
N. h/ps	2	2	4
INTERVENTION	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE
EXPECTED REPLACEMENTS			WIPER BRUSH ASPIRATOR + FILTERS DRAIN PUMP

6th YEAR			
	1st quarter	2nd quarter	3rd quarter
N. technicians	1	1	1
N. h/ps	2	2	5
INTERVENTION	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE
EXPECTED REPLACEMENTS			WIPER BRUSH ASPIRATOR + FILTERS FLANGE GASKETS DOOR GASKETS AIR/WATER TUBES- HANDPIECES

7th YEAR			
	1st quarter	2nd quarter	3rd quarter
N. technicians	1	1	1
N. h/ps	2	2	3
INTERVENTION	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE
EXPECTED REPLACEMENTS			WIPER BRUSH ASPIRATOR + FILTERS

8th YEAR			
	1st quarter	2nd quarter	3rd quarter
N. technicians	1	1	1
N. h/ps	2	2	4
INTERVENTION	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE
EXPECTED REPLACEMENTS			WIPER BRUSH ASPIRATOR + FILTERS FLANGE GASKETS DOOR GASKETS

9th YEAR			
	1st quarter	2nd quarter	3rd quarter
N. technicians	1	1	1
N. h/ps	2	2	4
INTERVENTION	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE
EXPECTED REPLACEMENTS			WIPER BRUSH ASPIRATOR + FILTERS AIR/WATER TUBES- HANDPIECES

	10th YEAR		
	1st quarter	2nd quarter	3rd quarter
<b>N. technicians</b>	1	1	1
<b>N. h/ps</b>	2	2	5
<b>INTERVENTION</b>	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE	PREVENTIVE MAINTENANCE
<b>EXPECTED REPLACEMENTS</b>			WIPER BRUSH ASPIRATOR + FILTERS FLANGE GASKETS DOOR GASKETS DRAIN PUMP